


# ZHANGJIAGANG MEDI MEDICAL EQUIPMENT CO., LTD.



## Specifications:

Dimension in total: 2180\*555\*750-1000mm

## Technical parameters:

1.HI-LOW: 750mm-1000mm

2.Head board fold upward $\geq 30^{\circ}$ and fold downward $\geq 90^{\circ}$

3.Back board fold upward $\geq 70^{\circ}$ and fold downward $\geq 15^{\circ}$

4.Leg board fold downward $\geq 90^{\circ}$ and fold outward $\geq 90^{\circ}$

5.Trendelenburg $\geq 25^{\circ}$ and reverse-trendelenburg $\geq 25^{\circ}$

6.Left and right tilting angle $\geq 20^{\circ}$

7.Waist board lifting $\geq 110$ mm

8.X-ray examination and photography function available, and C-shape arm applicable

9. Horizontal and longitudinal movement  $\geq 350$  mm

10. With controller inner constructed laterally for blocking of movements

## Structure: Stainless steel

This operation table can be for general surgical operation, thoracic operation, cardiovascular, extra-corporeal circulation, celiac, cholecystectomy nephrectomy, neck, cervical vertebra, lumbar spine, spine, back, urinary, gynecology, cranial nerve, surgery, rectum, and orthopedics lower limbs traction

Thoracic operation



Cranial nerve(Lying posture)



Cholecystectomy nephrectomy



Spine(Kneeling posture)



Spine(Leaned lying posture)



Cranial nerve(Sitting posture)



Cranial nerve(Leaned lying posture)



Upper extremity



Urinary and Gynecology



Hand



Rectum



Knee





Cardiovascular

Head and Neck



Eyes and Prosthetics