



张家港麦迪卡医疗器械有限公司
ZHANGJIANG MEDI MEDICAL EQUIPMENT CO., LTD.



CE ISO9001 ISO13485



Model: YA-XD1A ELECTRIC OPERATION TABLE

Size: (2150±50mm)*(500±20mm) *(750-1000)±50mm	
Application: For various surgical operations(neck, limbs etc.) also orthopedic operations with traction frame	
Functions:	
1.HI-LOW:(750-1000)±50mm	
2.Head board fold upward≥30°and fold downward≥90°	
3.Back board fold upward≥45°and fold downward≥15°	
4.Leg board fold downward≥90°and fold outward≥90°	
5.Trendelenburg≥25°and reverse-trendelenburg≥25°	
6.Left and right tilting angle≥15°	
7.Waist board lifting≥110mm	
8.X-ray examination and photography function available	
Structure: Stainless steel bottom base	
Accessories	
Imported electric pump for operation, Hand controller	
Power:AC220V 50HZ	Input power:300VA
Loading capacity:200kgs	
Packing: 148*76*95cm=1.1M3 G.W.:280kgs	